



QUICK START

- 1 Charge fully before first use.
- 2 Push the selected massage head firmly into the shaft.
- 3 Long-press the power button to turn on.
- 4 Short-press to choose a comfortable speed.
- 5 Move lightly over the muscle for about 60 seconds per area.
- 6 Long-press to turn off; wipe clean after use.

IMPORTANT

- Adult use only.
- Use only on dry, clean skin.
- Begin at a low speed with light pressure.
- Stop immediately if pain, numbness, tingling, dizziness, unusual discomfort, or skin irritation occurs.

PRODUCT SPECIFICATIONS

Model	MAS-G
Battery	1500 mAh
Rated power	25 W
Operating time	Approx. 2 hours
Charging time	Approx. 3 hours
Dimensions	130 x 155 x 55 mm
Product weight	0.4 kg
Housing material	ABS

OUTSIDE - FOLD-IN PANEL



CARE AND STORAGE

- Turn off and unplug before cleaning. Wipe with a soft, slightly damp cloth and dry thoroughly.
- Do not immerse in water or use abrasive cleaners or solvents.
- Store in a cool, dry place away from direct sunlight and excessive heat.
- A small amount of lubricant may be present near the shaft to reduce wear. Wipe away any excess before use.

LIMITED WARRANTY

This product includes a 12-month limited warranty from the original purchase date and is not transferable. Proof of purchase is required. The warranty excludes normal wear, accidental damage, misuse, improper transport, unauthorized repair or modification, and use contrary to these instructions.

SCAN FOR PRODUCT DETAILS



DISTRIBUTED BY UNITED YORAM INC.

Brooklyn, NY 11236

drfranklyns.com

support@drfranklyns.com | [@drfranklyns](https://twitter.com/drfranklyns)

OUTSIDE - BACK PANEL



DEEP TISSUE MASSAGE GUN

USER MANUAL

MODEL MAS-G



RECOVERY • RELAXATION • WELLNESS

READ INSTRUCTIONS COMPLETELY BEFORE USE

OUTSIDE - FRONT COVER

PRODUCT OVERVIEW



CONTROLS AND CHARGING

- Long-press the top power/control button to turn the device on or off.
- Short-press while operating to cycle through the available speed levels.
- The display indicates current operating information.
- Charge fully before first use. A first charge may take approximately 3 hours.
- Use a compatible 5V DC USB power source and the supplied charging cable.
- The device cannot be operated while charging. Disconnect before use.

MASSAGE HEADS



SPHERICAL HEAD

Large muscle groups: arms, shoulders, back, glutes, thighs, and calves.



U-SHAPED HEAD

Muscles beside the spine, calves, and Achilles area. Never apply directly over the spine.



BULLET HEAD

Focused pressure for palms, soles, and small muscle areas. Avoid joints and bones.



FLAT HEAD

Broad contact for general muscle relaxation and post-exercise massage.

GOOD TECHNIQUE

- Keep the head moving across the muscle.
- Use about 60 seconds per area, then reassess comfort.
- Check the treated area frequently; excessive pressure may cause bruising or soreness.

IMPORTANT SAFETY INFORMATION

- Keep away from children.
- Do not use on injured, swollen, inflamed, infected, numb, or painful areas.
- Never use on the head, throat, front or sides of the neck, bones, joints, spine, or internal-organ areas.
- Keep fingers, hair, loose clothing, and objects away from the moving shaft and ventilation openings.
- Never insert objects into the vents. Do not immerse the device in water or allow liquid to enter it.
- Do not drop, misuse, alter, repair, operate, or charge the device unattended.
- Consult a healthcare professional before use if pregnant; if you have a pacemaker or implanted device; take blood-thinning medication; have a circulatory or neurological condition; or recently had surgery.

BEFORE EVERY USE

- Inspect the device, cable, shaft, and attachment for damage.
- Insert the selected attachment firmly. Turn the device off before changing heads.
- Do not use if any component is damaged or operating abnormally.