

QUICK REFERENCE

Brand / Model	MagniPulse PEMF-01
Frequency	27.12 MHz
Automatic shutoff	After 15 minutes
Pulse rate	7 bursts per second
Power source	CR2025 button cell
Treatment area	110 cm²
Expected service life	Approximately 250 hours
Weight	Approximately 10 g

WHEN TO CONTACT A PROFESSIONAL

Stop use and seek medical advice if pain worsens, new symptoms appear, skin becomes irritated, or discomfort does not improve after 7 days.

DISPOSAL

This product contains electronic components and an integrated battery. Do not place it in fire or ordinary household trash where local e-waste collection is available. Follow local regulations.

CUSTOMER SUPPORT

Distributed by: _____
Website: _____
Email / Phone: _____

IMPORTANT INFORMATION

Read this manual completely before use. Keep it for future reference. This device is intended for single-user, non-sterile, external application. It is not a substitute for diagnosis or professional medical care.

IMPORTANT

Use only for the labeled intended purpose and according to the instructions supplied with the finished product.

SYMBOL GUIDE

- ! Caution - review warnings
- i Consult instructions for use
- BF Type BF applied part, if shown on product
- 250h Approximate total operating life
- LOT Batch / lot identification
- SN Serial number, if assigned
- Dispose through appropriate e-waste stream
- Keep dry

DOCUMENT CONTROL

Model: PEMF-01 | Revision: 1.2 | Issue date: September 2026

MagniPulse

PULSED ELECTROMAGNETIC FIELD THERAPY



MODEL PEMF-01

250-HOUR DEVICE

USER MANUAL

15-minute automatic shutoff • Portable • External use

1. ABOUT MAGNIPULSE

MagniPulse is a compact, battery-powered pulsed electromagnetic field (PEMF) device designed for convenient placement over an area of localized musculoskeletal discomfort. The treatment field is centered within the loop.

2. INTENDED USE

For adjunctive, temporary relief of minor localized musculoskeletal aches and pain. Use only on intact skin and only as described in this manual.

3. BEFORE FIRST USE

- Inspect the module, loop and covering. Do not use if cracked, cut, loose or wet.
- Confirm the skin is clean, dry and free of oils or lotion.
- Choose a wrap, bandage or medical-grade adhesive that holds the loop securely without restricting circulation.
- Record the date first activated. The device is powered by a CR2025 button cell.

4. CONTROLS AND INDICATORS



Button: press and hold approximately 1-2 seconds to switch on or off. Indicator: illuminated while operating. The unit switches off automatically after 15 minutes. Treatment area: 110 cm² within the loop.

5. HOW TO USE

- 1

POSITION
Center the loop directly over the area of discomfort. Keep the module and loop as close to the body as practical.
- 2

SECURE
Use a suitable wrap, bandage, garment or adhesive. The fit should be stable but comfortable; do not impair circulation.
- 3

TURN ON
Press and hold the button for approximately 1-2 seconds, then release. Confirm the indicator illuminates.
- 4

TREAT
A treatment cycle runs for 15 minutes and then switches off automatically. Recheck the skin and fit between cycles. Follow professional advice for repeat use.
- 5

TURN OFF
To stop early, press and hold the button for approximately 1-2 seconds. Confirm the indicator turns off before storing.

PLACEMENT TIPS

- Place the loop over, not merely near, the intended treatment area.
- Light clothing may be worn between the device and skin if the device remains securely positioned.
- Move the device only after switching it off. Inspect the skin before reapplying.

250-HOUR OPERATING LIFE
The CR2025 button cell provides an expected service life of approximately 250 hours. Actual life may vary with use and storage conditions.

6. DO NOT USE

- Directly over a pacemaker, implanted defibrillator, deep-brain stimulator, nerve stimulator or other active implanted electronic device.
- If pregnant or if you may be pregnant, unless specifically cleared by your healthcare professional.
- Over open wounds, infected, irritated or severely damaged skin.
- For sudden, severe or unexplained pain, chest or abdominal pain, cancer-related pain, or pain of unknown cause.
- On a child, or on anyone unable to understand the instructions or report discomfort, without professional supervision.

7. WARNINGS AND PRECAUTIONS

- Keep dry. Do not wear while bathing, showering or swimming.
- Do not crush, puncture, open, modify, heat or burn the device.
- Keep away from children and pets. Small parts may present a choking hazard.
- Do not share between users. The device is supplied non-sterile.
- Consult a healthcare professional before use if you are under medical care, have impaired sensation or circulation, or have concerns about electromagnetic therapy.

8. CLEANING AND STORAGE

Switch off. Wipe with a soft cloth lightly dampened with water and mild soap. Do not immerse, soak, sterilize or use solvents. Allow to dry fully. Store flat in a clean, dry place away from sunlight, excessive heat, moisture and sharp objects.

9. TROUBLESHOOTING

- Indicator does not turn on**
Repeat the 1-2 second press. Inspect for damage. If the device has reached its operating-life limit, discontinue use and dispose of it properly.
- Device will not stay in place**
Clean and dry the skin; replace the wrap or adhesive; do not overtighten.
- Skin irritation occurs**
Remove the device and adhesive. Wash the area and stop use until irritation resolves; seek advice if persistent.