

CARE AND PRODUCT LIFE

CLEANING

Follow the care label. Wipe surfaces with a soft, damp cloth and mild soap; air dry completely. Do not use bleach, strong solvents, or direct heat.

INSPECT BEFORE USE

Check the brace and closures for cracks, damage, loose parts, or loss of support. Stop using and replace if damaged.

EXPECTED LIFESPAN

Approximately three years under normal use and care. Actual life depends on wear, frequency of use, and condition. Replace sooner if needed.

CUSTOMER SUPPORT

ONE-YEAR LIMITED WARRANTY

This product carries a one-year limited warranty from the original purchase date. Contact customer support for warranty assistance and keep your proof of purchase.

CONTACT

Distributed by United Yoram Inc.
Brooklyn, NY 11236
drfranklyns.com
support@drfranklyns.com
Instagram: @drfranklyns

PRODUCT DETAILS

Scan the code below for more product information.



CERVICAL TRACTION NECK BRACE

SKU CC-03 | USER GUIDE



Read before use. Keep for reference.

PURPOSE AND USE

INTENDED PURPOSE

The cervical neck brace is intended to provide adjustable external support to the neck. It may be used as part of a healthcare professional's plan for neck discomfort or recovery.

BENEFITS

Designed to support neck positioning and help limit unwanted movement while worn.

BEFORE FITTING

For injury, whiplash, or after surgery, obtain professional fitting and wear instructions. Do not use as emergency stabilization for a suspected spinal injury.

FITTING AND WEAR

FIT WITH GUIDANCE

Have a healthcare professional demonstrate the correct placement, adjustment, and removal for your condition. Follow the clinician's prescribed position and wear schedule.

CHECK THE FIT

The brace should feel secure without pressing on the throat or restricting breathing or swallowing. Never force the neck into position or tighten to the point of pain.

DURING USE

Keep the head in the position advised by your clinician. Do not sleep, drive, or change prescribed settings unless your clinician says it is safe.

SAFETY INFORMATION

STOP AND SEEK ADVICE

Stop use and contact a healthcare professional if pain increases or you develop numbness, tingling, dizziness, headache, skin irritation, or unusual symptoms.

URGENT SYMPTOMS

Seek urgent medical attention for difficulty breathing or swallowing, new weakness, or severe symptoms.

IMPORTANT

Do not apply over an open wound unless directed. Do not alter the brace. If prescribed after surgery or injury, do not remove it contrary to your clinician's instructions.

TRACTION

Use any traction adjustment only as directed by a qualified healthcare professional. Do not force or over-adjust the brace.