

IMPORTANT SAFETY

- Use only as directed by a healthcare professional.
- Do not overtighten or block breathing, swallowing, or circulation.
- Remove immediately for increased pain, numbness, tingling, weakness, dizziness, breathing difficulty, or skin irritation.
- Do not drive, operate machinery, or sleep in the collar unless specifically cleared by a clinician.

BEFORE EACH USE

- Inspect the foam, fabric, seams, and closure.
- Make sure the collar is clean, dry, and free of damage.
- Confirm the prescribed wear schedule and fit instructions.

STOP USE

Remove the collar and contact a clinician if symptoms worsen, new neurological symptoms develop, or the collar causes pressure injury.

PRODUCT	Soft Cervical Collar
SKU	CC-1
SIZE	Adjustable: 15-20 in
MATERIAL	Latex-free foam / knit
WARRANTY	1 year
EXPECTED LIFE	3 years

INCLUDED

1 Adjustable Soft Cervical Collar

SCAN FOR PRODUCT DETAILS



Product information, support, and updates.

drfranklyns.com



Soft Cervical Collar Neck Brace

GENTLE SUPPORT AND STABILIZATION



USER MANUAL | CC-1

READ ALL INSTRUCTIONS BEFORE USE

HEALTH | SUPPORT | RECOVERY

PRODUCT OVERVIEW

The Dr. Franklyn's Soft Cervical Collar is designed to gently support and stabilize the cervical spine while limiting excessive neck movement. The contoured foam supports the chin and jaw, while the breathable cover and adjustable closure help provide a comfortable fit.



KEY FEATURES

- Soft, high-density latex-free foam support.
- Contoured profile supports the chin and jaw.
- Breathable knit fabric cover.
- Adjustable hook-and-loop closure.
- Lightweight design for prescribed daily use.

HOW TO FIT

1

PREPARE

Sit or stand upright with the head in a neutral, forward-facing position.

2

OPEN

Release the hook-and-loop strap and identify the higher contoured area for the chin.

3

POSITION

Center the collar under the chin. The lower edge should rest comfortably on the upper chest.

4

WRAP

Bring both ends around the neck without twisting the foam or fabric.

5

SECURE

Fasten the closure at the back or side as designed. Adjust until snug and supportive, never tight.

6

CHECK

You should breathe and swallow normally. Keep the chin centered and confirm the collar does not rub or pinch.



REAR CLOSURE

Overlap the ends evenly and secure the strap so the collar stays level. Recheck fit after movement.

WARNINGS AND PRECAUTIONS

ASK A CLINICIAN BEFORE USE IF YOU HAVE:

- A suspected fracture, dislocation, spinal instability, or recent neck trauma
- Recent cervical surgery or an open wound near the collar area
- Difficulty breathing or swallowing, severe osteoporosis, or neurological symptoms
- Skin sensitivity, impaired sensation, or circulation concerns

CLEANING AND STORAGE

- Fasten the closure before cleaning. Follow the sewn-in care label.
- When machine washing is permitted, use cold water, a gentle cycle, and a laundry bag.
- Air dry completely. Do not bleach, iron, or use high heat.
- Store clean, dry, and flat or loosely shaped away from sunlight and sharp objects.

SKIN AND FIT CHECKS

- Check the skin regularly, especially under the chin and along collar edges.
- Reposition and refasten if the collar shifts, folds, or causes pressure.
- Replace if the foam stays compressed, the fabric tears, or the closure no longer holds securely.

WEAR SCHEDULE

Wear only for the duration and activities directed by your healthcare professional. Do not extend wear time or use during sleep without approval.