

PRODUCT OVERVIEW

The Dr. Franklyn's CCW-1 combines intermittent pneumatic compression with a reusable cold-therapy gel pack. The universal wrap provides controlled inflation and release to simulate kneading and stroking of tissue, while localized cold therapy may help temporarily relieve minor muscle aches and support recovery.

FEATURES

- Five selectable pressure settings from 20 to 100 mmHg.
- Two air ports for use with one or two compatible wraps.
- Reusable non-toxic gel pack for localized cold therapy.
- Operates with two user-installed protected flat-top 18650 lithium-ion batteries; batteries are not included.
- Universal wrap for multiple treatment areas.

PACKAGE CONTENTS

- CCW-1 compression controller.
- Universal pneumatic wrap.
- Reusable cold-therapy gel pack.
- Air tubing, USB charging cable, and screwdriver.
- User manual. Batteries are not included.

QUICK REFERENCE

Model	CCW-1
Therapy	Intermittent compression + cold
Pressure range	20-100 mmHg
Pressure levels	5
Air ports	2
Power	Two protected flat-top 18650 Li-ion batteries (not included)

INTENDED USE AND MEDICAL PURPOSE

For adult home use to apply intermittent external compression, with optional cold therapy, to an intact lower leg or calf as directed. Provides temporary relief of minor muscle soreness, fatigue, discomfort, and localized swelling associated with activity or recovery. It is not intended to diagnose or treat a blood clot or vascular disease.

HOW IT WORKS

The controller cyclically inflates and releases the wrap at the selected pressure level. A reusable chilled gel pack may simultaneously provide localized cold therapy.

BEFORE FIRST USE

Read the entire manual. Install two approved protected flat-top 18650 batteries, confirm correct polarity, and inspect all components before operation.

1 Freeze the gel pack

Place the gel pack in a freezer for approximately 2 hours at -18°C ± 2°C (0°F ± 4°F).

2 Insert the gel pack

Lay the universal wrap flat and secure the frozen gel pack inside the designated wrap pocket.

3 Protect the skin

Always place one to two layers of clothing or a cloth barrier between the wrap and the skin.

OPERATING INSTRUCTIONS

- 1 Position and secure the wrap over the treatment area. Keep the fit snug, never excessively tight.
- 2 Connect the wrap tubing firmly to air port P1 or P2.
- 3 Press and hold the Power button to turn the controller on.
- 4 Press the pressure/gear control to choose a comfortable level from 20 to 100 mmHg. Begin low.
- 5 Press the Mode button to choose the desired intermittent compression program.
- 6 Monitor skin and comfort throughout the session.
- 7 After use, turn the controller off, loosen the wrap, and disconnect the tubing.

COLD-THERAPY LIMIT

Do not use the frozen gel pack for more than 30 minutes per session. Remove it before continuing compression-only therapy.

HOME AND REPEATED USE

Reusable by one person. Inspect the wrap, air bladder, tubing, ports, controller, and cold pack before every home session.

IMPORTANT

Do not use for suspected deep vein thrombosis, severe arterial disease, uncontrolled heart failure, active infection, open wounds, or unexplained calf pain or swelling without immediate medical evaluation.

WARRANTY AND EXPECTED LIFE

One-year limited warranty from the original purchase date. Expected product life: three years under normal use, inspection, care, cleaning, and storage. Stop using sooner if damaged, worn, contaminated, leaking, or unable to perform as intended. Retain proof of purchase.



MODEL CCW-1


COLD COMPRESSION THERAPY WRAP

Intermittent air compression + cold therapy

READ INSTRUCTIONS COMPLETELY BEFORE USE

WRAP PLACEMENT

The universal wrap may be positioned around the shoulder, waist/lower back, thigh, knee, calf, forearm, or elbow. Do not place it over the neck, chest, head, open wounds, or an area with impaired sensation.



PRESSURE & PORTS

- Start at the lowest comfortable pressure and increase only as needed.
- Use P1 for one wrap. P2 may support a second compatible wrap.
- Never kink, pinch, or sharply bend the air tubing.
- If numbness, throbbing, pulsation, pain, or unusual discomfort occurs, stop immediately and loosen the wrap.

IMPORTANT SAFETY INFORMATION

Consult a qualified healthcare professional before use if you have a medical condition or are unsure whether compression or cold therapy is appropriate for you.

DO NOT USE

- If deep vein thrombosis, thrombophlebitis, pulmonary edema, heart failure, severe arteriosclerosis, ischemic vascular disease, or an active infection is suspected or present.
- Over open wounds, recent skin grafts, dermatitis, gangrene, massive edema, or severe deformity.
- If you have neuropathy, extreme insensitivity to pain, significant vascular impairment, or a condition for which increased circulation or cold exposure is undesirable.
- On a limb where the wrap would interfere with vein ligation or another prescribed treatment.
- On children or anyone unable to communicate discomfort unless specifically directed by a healthcare professional.

WARNINGS & CAUTIONS

- Never place the frozen gel pack directly against exposed skin.
- Stop use if swelling, burning, excessive redness, numbness, skin discoloration, or increasing pain develops.
- Keep the controller, charger, and connections dry.
- Do not disassemble, modify, drop, crush, or expose the controller to extreme heat.
- Do not use near flammable anesthetic mixtures, oxygen-rich environments, or nitrous oxide.
- Use only compatible Dr. Franklyn's wraps, tubing, charger, and accessories.


CARE & CLEANING

- Turn off and disconnect the controller before cleaning.
- Wipe the controller with a soft, lightly damp cloth. Never immerse it.
- Remove the gel pack before cleaning the wrap. Wipe the wrap gently and allow it to air dry completely.
- Do not use bleach, abrasive cleaners, solvents, machine washing, or tumble drying unless the product label specifically permits it.

STORAGE

Store the controller, wrap, gel pack, and tubing in a clean, dry place away from direct sunlight, extreme heat, moisture, children, pets, and sharp objects. Do not store the tubing kinked or under heavy objects.

BATTERY INSTALLATION



- Disconnect the USB cable. Use the supplied screwdriver to loosen and remove the battery-cover screw.
- Slide the upper battery cover upward and remove it.
- Insert two matching protected flat-top 18650 rechargeable batteries, following the + and - polarity symbols.
- Reinstall the cover and securely tighten the screw before use.

BATTERIES REQUIRED - NOT INCLUDED

Protected flat-top 18650 | 1800-2600 mAh each | 65-68.5 x 18-18.6 mm

CUSTOMER SUPPORT

Distributed by United Yoram Inc. | Brooklyn, NY 11236
support@drfranklyns.com | drfranklyns.com
Instagram: @drfranklyns

DOCUMENT CONTROL

Brand: Dr. Franklyn's | Model: CCW-1 | Revision: 1.1 | September 2026

SCAN FOR PRODUCT DETAILS



drfranklyns.com

CCW-1 | REV 1.0 | SEPTEMBER 2026 | PAGE 2 OF 2