

CARE AND PRODUCT LIFE

CLEANING

Switch off and disconnect power before cleaning. Remove the gel pack. Wipe the wrap and controller with a soft, slightly damp cloth; dry fully. Do not immerse the controller or wash the powered assembly.

GEL PACK CARE

Wipe clean with mild soap and water. Do not puncture, microwave, heat, or use a leaking pack. Keep away from children and pets; gel is not for consumption.

STORAGE AND INSPECTION

Store clean and dry, away from sunlight and heat. Inspect seams, closures, gel pack, and controller before each use. Replace damaged parts.

EXPECTED LIFESPAN

Expected lifespan: 2-3 years, depending on frequency of use, care, storage, and normal wear and tear. Replace sooner if damaged, leaking, or no longer functioning correctly.

CUSTOMER SUPPORT

ONE-YEAR LIMITED WARRANTY


One-year limited warranty from the original purchase date. Retain proof of purchase and contact support for warranty assistance.

CONTACT

Distributed by United Yoram Inc.
Brooklyn, NY 11236
drfranklyns.com
support@drfranklyns.com
Instagram: @drfranklyns

PRODUCT DETAILS

Scan for more information about this product.
Quote SKU CCW-3 when contacting support.



COLD COMPRESSION THERAPY WRAP

Reusable cooling + adjustable support

SKU CCW-3 | USER GUIDE



Read all instructions before use.

Keep this guide for future reference.

PURPOSE AND BENEFITS

INTENDED USE

Reusable cold therapy wrap with adjustable compression support for temporary relief of muscle and joint discomfort. Cold may help ease soreness and swelling.

AREAS OF USE

Designed for use on the knee, elbow, arm, thigh, calf, and lower back. Use only where the wrap fits securely without restricting circulation.

INCLUDED COMPONENTS

Multi-area wrap, reusable gel pack, and adjustable straps. Includes an on/off controller.

PROFESSIONAL GUIDANCE

After surgery or a significant injury, obtain advice before use. Ask a healthcare professional before use with poor circulation, reduced sensation, diabetes, cold sensitivity, or a history of blood clots.

PREPARATION AND APPLICATION

1 | PREPARE THE GEL PACK

Follow the cooling instructions printed on the gel pack. Do not chill or freeze the controller or powered wrap. Do not use the pack if damaged or leaking.

2 | POSITION THE PACK

Place the cooled pack in the wrap's designated pocket or holder. Keep a fabric barrier between the cold pack and your skin. Never apply a frozen pack directly to bare skin.

3 | SECURE THE WRAP

Place over the intended area. Fasten straps evenly for a comfortable fit. Do not pull tight enough to cause numbness, tingling, or discoloration.

4 | POWERED OPERATION

Use the marked on/off control only after fitting. Follow the markings supplied with your unit. Ask support for instructions if controls or power requirements are unclear. Use only the power supply specified for your unit.

SAFE TREATMENT AND REMOVAL

DURING TREATMENT

Remain awake and able to monitor comfort. Check the area regularly. Follow your clinician's recommended session length; stop sooner if uncomfortable. Allow normal skin temperature and sensation to return between sessions.

STOP USE IMMEDIATELY

Switch off and loosen or remove the wrap for pain, burning, excessive cold, numbness, tingling, unusual swelling, or pale, blue, or blotchy skin. Seek medical advice if symptoms persist.

DO NOT USE

Do not use while sleeping, over open wounds, or on areas with reduced sensation without professional guidance. Do not use for a suspected blood clot or an unexplained hot, swollen limb; seek medical assessment.

AFTER USE

Turn off before loosening straps. Remove the gel pack, inspect for damage, and let the wrap dry before storage. Keep electrical components away from water.