

IMPORTANT SAFETY

- Use only as directed by a qualified healthcare professional.
- Start with gentle traction and a short session.
- Stop immediately for pain, dizziness, numbness, tingling, weakness, nausea, headache, or visual changes.
- Never sleep while using the device. Keep away from children.

BEFORE EACH USE

- Inspect the frame, cushion, strap, tubing, pump, and valve.
- Do not use damaged, leaking, loose, or worn parts.
- Use only on a firm, stable, flat surface.

STOP USE

Release all air immediately if discomfort or neurological symptoms occur. Seek medical advice if symptoms persist.



PRODUCT INFORMATION

PRODUCT	Cervical Traction Device
SKU	CT-1
WARRANTY	1 year
EXPECTED LIFE	3 years

CARE AND STORAGE

Wipe with a soft cloth lightly dampened with mild soap and water. Do not immerse. Dry completely and store deflated in a cool, dry place away from heat, sunlight, and sharp objects.

SCAN FOR PRODUCT DETAILS



Instructions, product information, and support.

drfranklyns.com

United Yoram Inc. | Brooklyn, NY 11236 | support@drfranklyns.com



Cervical Traction Device

GENTLE ADJUSTABLE NECK TRACTION




[USER MANUAL | CT-1](#)

READ ALL INSTRUCTIONS BEFORE USE

HEALTH | SUPPORT | RECOVERY

PRODUCT OVERVIEW

The Dr. Franklyn's Cervical Traction Device is a manually inflated, non-powered device designed to apply gentle, adjustable cervical traction while lying down. Use only as recommended by a qualified healthcare professional.




DEVICE COMPONENTS

- Forehead securing strap
- Inflatable cervical traction cushion
- Traction frame
- Hand inflation pump
- Air-release valve and tubing

HOW TO USE

- 1 PREPARE**
Place on a firm bed or therapy table. Loosen the forehead strap and mostly deflate the cushion.
- 2 POSITION**
Lie on your back. Center the back of your neck on the cushion with shoulders relaxed and head aligned.
- 3 SECURE**
Place the strap comfortably across the forehead. Fasten securely, but not tightly. Keep it away from the eyes and throat.
- 4 INFLATE SLOWLY**
Close the release valve. Squeeze the hand pump gradually until you feel a gentle, comfortable stretch.
- 5 BEGIN CONSERVATIVELY**
Follow the duration and traction level recommended by your healthcare professional. Begin with low pressure and a short session.
- 6 RELEASE AND EXIT**
Open the valve slowly. Fully deflate the cushion, unfasten the strap, and sit up gradually.



CORRECT POSITION

Rest the neck over the cushion while lying flat. Keep the head centered and shoulders relaxed.


WARNINGS AND PRECAUTIONS

ASK A CLINICIAN BEFORE USE IF YOU HAVE:

- Recent fracture, surgery, dislocation, or acute neck injury
- Spinal instability, severe osteoporosis, or spinal cord compression
- Rheumatoid arthritis of the neck, infection, malignancy, or inflammatory cervical disease
- Severe vascular disease, uncontrolled high blood pressure, or another condition that makes traction unsafe


IMPORTANT PRECAUTIONS

- Never exceed a comfortable traction level. More pressure is not necessarily more effective.
- Do not twist, turn, or lift the head while inflated. Breathe normally.
- Do not modify the device or use it for another purpose.
- Seek medical care for persistent or worsening symptoms.



STOP IMMEDIATELY

Release all air if pain, dizziness, numbness, tingling, weakness, nausea, headache, visual changes, or unusual symptoms occur.



SUPPORT

support@drfranklyn.com | drfranklyn.com