




CUSTOMER SUPPORT

Distributed by United Yoram Inc. | Brooklyn, NY 11236
drfranklyns.com
support@drfranklyns.com
Instagram: @drfranklyns

PRODUCT INFORMATION

Product: Dr. Franklyn's Portable Air Compression - Multi-Use
Model / SKU: DVT-1
Manual revision: 1.1 | September 2026
Limited warranty: 1 year
Expected product lifespan: 3 years

SCAN FOR PRODUCT DETAILS



drfranklyns.com

KEEP THIS GUIDE

Keep this manual for future reference. Follow the prescribing clinician's directions whenever they differ from general guidance in this booklet.

IMPORTANT SAFETY INFORMATION

The Dr. Franklyn's Portable Air Compression - Multi-Use provides intermittent air compression for wellness and recovery. It is not a substitute for medical evaluation or prescribed DVT prevention or treatment.

- Use only under clinician direction if you have or may have a blood clot, deep vein thrombosis, pulmonary embolism, severe vascular disease, or unexplained pain or swelling in an arm or leg.
- Do not use while the device is charging, while driving, or while sleeping.
- Do not over-tighten, puncture the air bladder, place it around the neck, chest, abdomen, or head, use near water, or expose it to flammable materials, excessive heat, or open flame.
- If a heating function is included and activated, avoid prolonged direct heat contact and check the skin frequently to prevent burns.
- Stop use if the cord is damaged, liquid enters the device, the device is dropped or damaged, or operation becomes abnormal.
- Do not disassemble, repair, or modify the controller, air bladder, battery, or charging components.

DO NOT USE / ASK A CLINICIAN

- Do not use during pregnancy or the postpartum period without medical approval.
- Ask a clinician before use with cancer, serious heart disease, abnormal blood pressure, serious liver disease, bone disease including Paget disease, skin allergy, or fragile/damaged skin.
- Ask a clinician before use if you have an implanted or personal medical device or cannot feel or communicate pressure, heat, or discomfort normally.

STOP USE

Stop immediately and seek medical advice for chest pain, shortness of breath, sudden limb swelling, severe or increasing arm or leg pain, numbness, skin discoloration, dizziness, burning, or any unusual reaction.



PORTABLE AIR COMPRESSION MULTI-USE

INTERMITTENT COMPRESSION FOR WELLNESS &
RECOVERY



MODEL DVT-1

USER MANUAL

READ ALL INSTRUCTIONS BEFORE USE

HEALTH | SUPPORT | RECOVERY

INTENDED USE

The Dr. Franklyn's Portable Air Compression - Multi-Use is intended to apply cyclic air pressure around an arm, leg, calf, or thigh to support circulation, relaxation, and recovery. Follow a healthcare professional's instructions when the device is part of a medical care plan.

IMPORTANT

This consumer device does not diagnose, treat, cure, or prevent deep vein thrombosis. Do not delay medical care for possible clot symptoms.

PRODUCT FEATURES

- Portable multi-use wrap for the arm, leg, calf, or thigh.
- Adjustable hook-and-loop straps for a secure, personalized fit.
- Integrated rechargeable lithium-ion battery.
- Selectable operating modes and automatic 15-minute session timer.
- USB charging with charging-status indicator.

TECHNICAL SPECIFICATIONS

Item	Specification
Product	Portable air compression - multi-use
Brand model	DVT-1
Rated input	5 V DC
Rated power	7.5 W
Default session	15 minutes
Power source	Built-in rechargeable Li-ion battery

BEFORE FIRST USE

- Inspect the wrap, air bladder, controller, charging cable, seams, and closures. Do not use damaged components.
- Charge fully using the supplied USB cable and a suitable 5 V USB power source. Do not operate while charging.
- The indicator flashes red while charging and turns off when charging is complete.

HOW TO OPERATE

1. Sit comfortably and relax the arm or leg being treated. Open the hook-and-loop straps.
2. Position the device around the selected arm, leg, calf, or thigh area with the control panel facing outward and the material lying smooth against the skin or a thin garment.
3. Keep the wrap away from the elbow and knee joint. Fasten the straps evenly so the device is secure without restricting circulation.
4. If supplied separately, connect the controller and air-bladder connection securely.
5. Short-press ON/OFF to begin. Use ON/OFF/MODE controls to select the available compression setting.
6. The device stops automatically after approximately 15 minutes. To stop earlier, press and hold ON/OFF for about 3 seconds.

DURING A SESSION

Compression should feel firm and rhythmic, not painful. Remain awake and regularly check the fingers or toes beyond the wrapped area for normal color, warmth, and sensation. Stop immediately for discomfort or abnormal symptoms.

CLEANING AND CARE

- Power off and disconnect the charging cable before cleaning.
- Do not wash, immerse, rinse, or soak the controller or assembled device.
- Wipe surfaces with a soft cloth lightly dampened with water or 70–75% alcohol, then allow to dry completely.
- Keep liquids out of the charging port, controller, tubing, and air connections.

TROUBLESHOOTING

Issue	Action
Weak compression	Reposition and refasten the wrap; select another available level.
Low sound or dim light	Battery may be low; recharge fully.
Will not turn on	Recharge; confirm the cable is disconnected before use.
Abnormal operation	Stop use, power off, and contact customer support.

STORAGE AND BATTERY

- Store in a cool, dry place away from children, moisture, direct sunlight, sharp objects, and excessive heat.
- For long-term storage, recharge approximately every two months.
- The device contains a lithium-ion battery. Do not place it in fire or household trash; follow local battery/electronic-waste rules.

WHAT IS INCLUDED

Portable air compression device with controller, supplied USB charging cable, and this user manual. Included accessories may vary by package configuration.