



Round Cervical Support Roll

NECK, BACK & POSITIONING COMFORT



INTENDED USE & PURPOSE

A cylindrical foam bolster intended to support comfortable positioning of the neck, lower back, knees or legs. It fills the space beneath a body curve to provide cushioning and help maintain a comfortable resting posture. Use for everyday comfort or clinician-directed positioning.

PRODUCT DETAILS

Non-powered foam bolster; round profile; removable, washable cover. Includes one cervical support roll with cover.

HOW TO USE

- 1. Neck:** Lie on your back and place the roll beneath the natural curve of your neck. Support your head as needed; avoid forcing the chin upward or toward the chest.
- 2. Lower back:** Place the roll horizontally between your lower back and the chair back. Sit back gently and adjust the height for comfort.
- 3. Knees or legs:** Place beneath the knees while lying on your back, or beneath the legs as directed. Keep pressure comfortable and evenly distributed.

IMPORTANT SAFETY

Use on a stable surface. Stop if pain, tingling, numbness or pressure increases. Seek professional advice after a neck injury or surgery, or if a specific position has been prescribed. This cushion does not immobilize the neck and should not replace a prescribed brace. Keep away from infants and unsupervised children.

CLEANING & CARE

Remove the cover and wash according to its care label. Do not machine wash or soak the foam. Spot clean the foam gently only as needed and let it dry fully. Store dry and away from direct heat. Replace if the foam or cover deteriorates.

PRODUCT SUPPORT

Distributed by United Yoram Inc. | Brooklyn, NY 11236
drfranklyns.com | support@drfranklyns.com
Instagram: @drfranklyns

