

EMS/TENS Compression Belt

ELECTRODE RETENTION & LIGHT WAIST SUPPORT



Belt shown with illustrative stimulation accessories.

INTENDED USE & PURPOSE

An adjustable belt accessory intended to help hold compatible EMS/TENS electrode pads in place around the waist or lower-back region while providing light compression. The belt itself does not generate electrical stimulation. Therapeutic use, settings and electrode placement depend on the separate stimulator.

PRODUCT DETAILS

Soft, breathable, flexible belt; adjustable fastening system. Includes one belt. Confirm compatibility with your stimulator and electrodes. Stimulator, leads and pads shown are illustrative and are not listed as included with EMS-B.

HOW TO USE

1. Read the separate stimulator and electrode instructions. Obtain advice on suitable electrode placement and treatment settings.
2. With the stimulator switched off, apply its electrodes as instructed to clean, intact skin. Connect the leads as directed.
3. Wrap the belt gently around the waist over the electrode area. Fasten snugly without restricting breathing, pinching cables or compressing the controller.
4. Operate the stimulator only as directed by its manufacturer. Switch it off before adjusting pads or removing the belt.

IMPORTANT SAFETY

Follow every warning in the stimulator manual. Do not use electrical stimulation with a pacemaker or other implanted electronic device unless specifically cleared by the treating specialist. Seek medical advice before use during pregnancy or with epilepsy. Never place electrodes across the chest, on the front of the neck, head, broken skin or numb areas. Do not use stimulation while driving, asleep, bathing or operating machinery. Stop for burning, painful stimulation or skin irritation.

CLEANING & CARE

Switch off and remove the stimulator, leads and electrodes before cleaning the belt. Follow the belt care label; hand clean gently with mild soap if permitted and air dry fully. Do not wash electrical parts or immerse them. Inspect the belt and fasteners before reuse.

PRODUCT SUPPORT

Distributed by United Yoram Inc. | Brooklyn, NY 11236
drfranklyns.com | support@drfranklyns.com
Instagram: @drfranklyns

