



Adjustable Hip & Thigh Support Brace

SOFT HIP, GROIN & UPPER-THIGH SUPPORT



Support shown fitted to the hip and upper thigh.

INTENDED USE & PURPOSE

A soft external support intended to provide gentle compression and support around the hip, groin and upper thigh. May be used for clinician-directed support during recovery from a soft-tissue strain or injury. After surgery, use only if approved by the treating clinician.

PRODUCT DETAILS

Flexible fabric wrap with adjustable waist and thigh hook-and-loop straps; lightweight, breathable design. Non-powered. Includes one hip-and-thigh support brace. This soft wrap does not lock a prescribed hip-abduction angle.

HOW TO USE

1. Open the waist and thigh straps. Position the connecting side panel along the outside of the affected hip and upper thigh.
2. Wrap and fasten the waist section securely above the hip. Smooth the material so there are no folds or pressure points.
3. Wrap the thigh section around the upper thigh and fasten the straps evenly. Avoid bunching in the groin.
4. Check comfort while sitting and standing. The fit should be snug without restricting circulation. Follow your clinician's wear and activity instructions.

IMPORTANT SAFETY

Do not use this wrap as a substitute for a prescribed rigid hip orthosis or dislocation precautions. Obtain medical evaluation after trauma, suspected fracture or dislocation. Do not apply over open wounds or an unprotected incision. Loosen and stop for numbness, swelling, discoloration, cold skin or increased pain. Seek advice if symptoms persist. Inspect skin regularly.

CLEANING & CARE

Close the hook-and-loop fasteners before cleaning. Follow the care label. Hand wash gently in cool water with mild soap if permitted; air dry fully. Do not bleach, iron or tumble dry. Replace if fabric, straps or fasteners no longer provide a secure fit.

PRODUCT SUPPORT

Distributed by United Yoram Inc. | Brooklyn, NY 11236
drfranklyns.com | support@drfranklyns.com
Instagram: @drfranklyns

