



Ergonomic Memory Foam Seat Cushion

SEATED COMFORT & TAILBONE PRESSURE RELIEF



INTENDED USE & PURPOSE

A contoured seat cushion intended to distribute sitting pressure and improve everyday seated comfort. Its rear U-shaped cutout reduces direct contact beneath the tailbone. Suitable for compatible office, home and vehicle seating, or wheelchair seating when approved by the seating clinician.

PRODUCT DETAILS

High-density memory foam; U-shaped rear cutout; breathable, removable cover; portable, non-powered construction. Includes one seat cushion and cover. SKU: SC-01.

HOW TO USE

1. Place the cushion flat on a firm, stable seat. Position the U-shaped cutout at the rear, beneath the tailbone.
2. Sit centrally with your hips supported on both sides of the cutout. Keep your feet supported and adjust the chair height if necessary.
3. Check that the cushion does not slide or interfere with the chair, seat belt or vehicle controls.
4. Change position and take regular movement breaks. Recheck cushion placement after standing or transferring.

IMPORTANT SAFETY

This cushion is not a substitute for a prescribed pressure-injury prevention system. If you have reduced sensation, fragile skin, a pressure injury or limited ability to reposition, ask your clinician to assess suitability. Stop for skin redness, discomfort or numbness. In a vehicle, confirm safe posture and control reach before driving.

CLEANING & CARE

Remove the cover and follow its care label. Do not soak or machine wash the foam. Spot clean the foam gently as needed and dry fully before refitting the cover. Store dry, away from heat and heavy loads.

PRODUCT SUPPORT

Distributed by United Yoram Inc. | Brooklyn, NY 11236
drfranklyns.com | support@drfranklyns.com
Instagram: @drfranklyns

