



Lumbar Support Cushion

LOWER-BACK SUPPORT & SEATED COMFORT



INTENDED USE & PURPOSE

A contoured foam back cushion intended to support the natural inward curve of the lower back during sitting. It helps fill the gap between the lumbar region and the chair back to encourage a comfortable seated posture and reduce unsupported sitting strain. Use for daily comfort or clinician-directed positioning.

PRODUCT DETAILS

High-density supportive foam; shaped lumbar profile; breathable fabric cover; adjustable securing strap; non-powered construction. Includes one lumbar support cushion with strap. For compatible office, home, vehicle or wheelchair seating.

HOW TO USE

1. Place the cushion upright against the chair back, with the curved support surface facing your lower back.
2. Position the cushion at the natural inward curve of your lumbar region. Wrap the securing strap around the chair back and adjust it so the cushion stays in place.
3. Sit back gently with your hips toward the back of the seat and your feet supported. Adjust the cushion height until it supports you without pushing you forward excessively.
4. Recheck comfort and placement regularly. Change position and take movement breaks. Adjust a vehicle cushion only while parked.

IMPORTANT SAFETY

Use on a stable, compatible seat. Do not obstruct a seat belt, airbag, wheelchair support or seat control. In a vehicle, confirm safe posture, visibility and pedal reach before driving. Stop if pain, numbness, tingling or skin irritation increases. Ask a healthcare professional about use after injury or surgery. This cushion does not immobilize the spine or replace a prescribed brace.

CLEANING & CARE

Follow the cover care label. Wipe with a soft, lightly damp cloth and allow to dry fully. If removing the cover for cleaning, follow its label instructions. Do not soak or machine wash the foam. Store dry, away from direct heat and heavy loads. Replace if the foam loses support or the strap or cover is damaged.

PRODUCT SUPPORT

Distributed by United Yoram Inc. | Brooklyn, NY 11236
drfranklyns.com | support@drfranklyns.com
Instagram: @drfranklyns

