



Premium Memory Foam Lumbar Support Cushion

CONTOURED SUPPORT FOR SEATED COMFORT



INTENDED USE & PURPOSE

An external seating support intended to cushion the lower back and support its natural curve during sitting. Slow-rebound memory foam contours to the body to improve seated comfort at home, at a desk or in a compatible vehicle seat.

PRODUCT DETAILS

High-density, slow-rebound memory foam; contoured lumbar profile; breathable perforated cover; adjustable securing strap. Non-powered. Includes one lumbar support cushion.

HOW TO USE

1. Place the cushion upright against the chair back, with its contoured surface facing you.
2. Position the fullest support at the natural inward curve of your lower back. Adjust the height before securing the strap around the chair back.
3. Sit with your hips toward the back of the seat and feet supported. The cushion should fill the gap without pushing you forward excessively.
4. Recheck comfort and position during use. Take regular movement breaks. In a vehicle, adjust the cushion only while parked.

IMPORTANT SAFETY

Use only on stable, compatible seating. Do not obstruct a seat belt, airbag or seat control. Remove the cushion if it prevents safe driving posture or pedal reach. Stop for increased pain, numbness or skin irritation. Obtain positioning advice after surgery or for persistent back symptoms. This is a cushion, not a spinal immobilizer.

CLEANING & CARE

Follow the cover care label. Wipe with a soft, lightly damp cloth and dry fully. Do not soak or machine wash the memory foam. Keep away from direct heat and store without heavy objects compressing the cushion.

PRODUCT SUPPORT

Distributed by United Yoram Inc. | Brooklyn, NY 11236
drfranklyns.com | support@drfranklyns.com
Instagram: @drfranklyns

