

PRODUCT INFORMATION

Keep this guide for future reference. Read all instructions and follow your healthcare professional's directions when the brace is prescribed.

CUSTOMER SUPPORT

Distributed by United Yoram Inc. | Brooklyn, NY 11236
drfranklyns.com | support@drfranklyns.com
Instagram: @drfranklyns

DOCUMENT CONTROL

Brand: Dr. Franklyn's | Model: LSO-2 | Revision: 1.0 | Issue date: September 2026

SCAN FOR PRODUCT DETAILS



drfranklyns.com

INTENDED USE AND MEDICAL PURPOSE

The LSO-2 is a reusable lumbar-sacral support intended to provide external compression and stabilization to the lower back. It may help manage lower-back strain, muscle fatigue, sciatica-related discomfort, disc-related symptoms, and clinician-directed rehabilitation or posture support.

SUPPORT PROVIDED

- Flexible stabilizing stays help support lumbar alignment while allowing functional movement.
- Dual compression straps allow the user to adjust support gradually and evenly.
- Breathable mesh helps reduce heat buildup during daily wear.
- Low-profile construction may be worn over or under clothing when comfortable.

DO NOT USE OR SEEK ADVICE

- Do not place over open wounds, infected skin, or an unprotected surgical incision.
- Seek professional advice after surgery, for suspected fracture, severe osteoporosis, major deformity, or persistent unexplained pain.
- Do not overtighten or continue use if symptoms increase or circulation or sensation changes.



HEALTH | SUPPORT | RECOVERY



MODEL LSO-2

ADJUSTABLE LSO
BACK BRACE

LOWER-BACK SUPPORT AND ADJUSTABLE COMPRESSION

READ INSTRUCTIONS COMPLETELY BEFORE USE

BEFORE FITTING

- Inspect the brace and closures. Ensure the stays are seated and the fabric is dry.
- Loosen the compression straps and position the brace over a thin garment when possible.

FITTING INSTRUCTIONS

- 1

CENTER

Place the rear support directly over the lower back with the bottom edge above the hips.
- 2

WRAP

Bring both main belt ends around the waist without twisting.
- 3

FASTEN

Secure the front closure so the brace is snug but comfortable.

FIT CHECK

The brace should be centered, level, and comfortably snug. It should support the lumbar region without restricting normal breathing or creating numbness or pressure points.

COMPLETE THE FITTING

- 4

TENSION

Pull both secondary straps forward with equal force.
- 5

ATTACH

Fasten the straps symmetrically to create comfortable compression.
- 6

RECHECK

Breathe and move normally. Reposition if the brace rolls, shifts, or pinches.

USE AND WEAR

- Follow a clinician's schedule when prescribed. Otherwise begin with shorter wear periods and check the skin.
- Do not sleep in the brace or use it during activities that worsen symptoms unless directed.

Remove the brace and seek medical advice for increasing pain, numbness, tingling, unusual swelling, skin breakdown, or pale, blue, cold, or discolored skin.

IMPORTANT SAFETY INFORMATION

- Do not rely on the brace to lift beyond your ability or to continue painful activity.
- Use under supervision if you have impaired circulation, sensation, skin integrity, or balance.
- Single-patient use only. Do not alter the stays or sew through structural panels.

CLEANING AND STORAGE

- Fasten all hook-and-loop closures before cleaning. Remove detachable rigid stays only when the care label permits.
- Hand wash in cool water with mild soap. Do not bleach, wring, machine dry, iron, or dry clean.
- Air dry completely away from direct heat and sunlight. Store flat or loosely folded in a clean, dry place.

WARRANTY AND EXPECTED LIFE

One-year limited warranty from the original purchase date. Expected product life: three years under normal fitting, inspection, care, cleaning, and storage. Replace sooner if damaged, stretched, or unable to provide the intended support. Retain proof of purchase.