

READ BEFORE USE

Safety first

Protect your eyes

Never look into the treatment head or aim it at anyone's eyes. Avoid mirrors and reflective surfaces. Use only manufacturer-specified eye protection matched to the device's wavelengths and output; sunglasses are not a substitute.

Get advice before treatment

Consult a healthcare professional if you are pregnant, have cancer or a suspicious lesion in the area, have reduced sensation, or take medicines that increase light sensitivity.

Use with care

Keep away from children. Do not use while sleeping, near water, or on damaged skin. Do not use a damaged device or modify the housing.

STOP IF DISCOMFORT OCCURS

Stop use for burning, irritation, unusual warmth, or increased pain. Seek medical advice for persistent symptoms. Seek prompt medical attention for any suspected eye exposure or vision change.

CARE & SUPPORT

Keep it ready

After each use

Switch off and disconnect the charging cable before cleaning. Wipe the housing with a soft, dry cloth. For treatment-window cleaning, contact the distributor for care instructions. Keep liquid out of openings.

Storage

Store dry and protected from heat, direct sunlight, moisture, and impact. Do not immerse, spray, or sterilize. Do not open the housing or attempt battery repairs.

Need help?

Stop use if the device, cable, or charging port is damaged, becomes unusually hot, or gives off an odor. Contact the distributor with SKU LT-2, your purchase details, and a description of the issue.

DISTRIBUTED BY

RAM Basics Inc.
Lakewood NJ 08701
Sales@rambasics.com

drfranklyns.com



SCAN FOR MORE
PRODUCT DETAILS



LASER THERAPY

Targeted muscle & joint comfort



USER MANUAL | SKU LT-2

Read before use. Keep for future reference.

What it is used for

Dr. Franklyn's Laser Therapy

Dr. Franklyn's Laser Therapy is intended to help manage muscle and joint discomfort, soreness, and stiffness. Treatment areas include the back, shoulders, knees, elbows, hips, hands, and neck, subject to the safety precautions in this manual and your clinician's guidance.

Role in a treatment plan

A clinician may consider light therapy as an additional treatment for pain or stiffness that limits movement or daily activity. Goals may include symptom relief and improved participation in exercise or rehabilitation. Results vary by condition, device, and treatment dose.

Medical necessity

Medical necessity is determined for the individual patient by the treating clinician. Documentation should explain the diagnosis, affected area, pain and functional limitations, treatments already tried, and why this device is appropriate within the care plan.

CLINICIAN TREATMENT PLAN

Record the intended goals, settings, treatment time, frequency, duration of care, safety precautions, and reassessment plan. Continued use should be based on the patient's response.

Specifications

Light wavelengths

650 nm red + 850 nm near-infrared

Light array

16 LEDs

Intensity adjustment

7 power levels

Timer

5 settings spanning 5-30 minutes

Stated output

9 W

Design

Handheld device with digital display

Package contents

Device, USB charging cable, user manual

Controls shown

Power: intensity selection.
Timer: session-time selection.
On/off: device operation.

Using the device

1 Inspect & charge

Check the housing, treatment window, cable, and port for damage. Charge on a dry, stable surface using the supplied cable and a source matching the device label. Disconnect the charging cable before treatment.

2 Prepare & protect

Clean and expose the treatment area as directed. Avoid mirrors and bystanders. Use the specified eye protection before emission. Do not look into the treatment head, even if some emitters appear unlit.

3 Select your treatment

Use only the intensity, time, frequency, and placement specified by your clinician for this device. The 5-30 minute timer range is a device feature, not a recommendation to treat for 30 minutes.

4 Treat, then switch off

Use the treatment distance or contact position specified by your clinician. Activate the device only when it is positioned away from eyes. Stay awake and attentive. Stop for discomfort. Stop emission before repositioning; switch off after the session.