



1 CARE AND STORAGE

- Switch off and disconnect power before cleaning.
- Wipe with a soft, slightly damp cloth. Never immerse the product or controls.
- Store clean and dry, away from direct sunlight, heat, moisture, and children.
- Stop using the product if its housing, cable, adapter, or treatment surface is damaged.

2 SUPPORT AND WARRANTY CARE

For assistance, provide model OPS-1, proof of purchase, and a description of the issue. Do not open or repair the product yourself.

3 BEFORE FIRST USE

- Place the pad on a stable chair or vehicle seat and secure all straps.
- Route cables away from pedals, seat tracks, and other driving controls.
- Use only the correct adapter and verify its voltage before connecting.

4 COMFORT FEATURES

- Targeted massage zones for neck, back, waist, hips, and legs.
- Multiple massage modes and strength settings.
- Optional gentle heat and magnetic massage features.
- Automatic timer and overheat protection.

IMPORTANT

Read all warnings before use. Begin with the lowest comfortable setting and stop if discomfort occurs.



MASSAGER SEAT

MODEL OPS-1

USER MANUAL | READ BEFORE USE



SCAN FOR PRODUCT DETAILS

Open the Dr. Franklyn's product page for current product information and support.

MODEL OPS-1

1 INTENDED USE

The Dr. Franklyn's OPS-1 Massager Seat is a seat-mounted vibration massager intended to provide temporary relaxation and soothing comfort for the neck, back, waist, hips, and legs. This product is not a child restraint or safety seat.

2 SAFE PLACEMENT

- Center the pad with the back section upright.
- Secure every strap before sitting.
- Confirm the pad cannot slide and no buckle, strap, or cable is beneath the user.
- The driver must not adjust the controller while the vehicle is moving.

RELAXATION AT HOME OR ON THE GO

Use the correct adapter for home, office, or vehicle operation. Always stop the vehicle safely before changing any setting.

3 OPERATING INSTRUCTIONS

- 1 Connect the controller and correct power adapter.
- 2 Press POWER to begin in the default massage mode.
- 3 Press MODE to cycle through massage patterns.
- 4 Press INTENSITY/STRONG to adjust strength.
- 5 Use ZONE controls, if present, for targeted areas.
- 6 Press HEAT to turn gentle warmth on or off.
- 7 Press POWER to stop, then disconnect the adapter.

4 AUTOMATIC OPERATION

Massage patterns may vary automatically. The unit shuts off after its programmed cycle. Allow it to cool before restarting.

COMFORT TIP

Begin at the lowest intensity. Sit upright against the cushion and increase the setting only when the massage feels comfortable.

! WARNINGS

- Do not use while sleeping, on wet skin, or in a bathroom. Keep the controller and connectors dry.
- Do not fold, crush, puncture, or place heavy objects on the pad. Never insert pins.
- Do not use on areas with reduced sensation, inflammation, injury, or impaired circulation.
- Ask a clinician before use if pregnant, using a pacemaker, recovering from surgery, or managing a serious condition.
- Stop immediately for pain, dizziness, numbness, burning, or excessive heat.

5 TROUBLESHOOTING

- No power: check the outlet, adapter, connectors, and vehicle accessory socket.
- Stops during use: the timer or overheat protection may have activated. Disconnect and cool.
- Controls do not respond: switch off, reconnect all plugs firmly, and restart.
- Unusual odor, noise, or damage: stop use and contact support.

WHEN TO STOP

Turn the unit off immediately if you feel pain, dizziness, numbness, burning, excessive heat, or unusual discomfort.