

Medical purpose & specifications

Dr. Franklyn's Deluxe Shoulder Immobilizer Brace with Adjustable Abduction Wedge (SS-04) is an external shoulder orthosis intended to support the upper extremity and restrict unwanted shoulder movement during clinician-directed recovery from injury or surgery.

CLINICAL PURPOSE & SUPPORT

For prescribed post-operative shoulder immobilization, including rotator cuff repair, and support after shoulder injury or dislocation. The forearm platform supports the elbow, forearm and wrist; the torso strap secures the assembly. The wedge maintains the prescribed arm position away from the torso.

PRODUCT SPECIFICATIONS

Model: SS-04

Use: Single patient; left or right arm

Materials: Foam padding, nylon, polyester and hook-and-loop fasteners

Components: Arm support platform, body panel, shoulder and waist straps, arm straps and adjustable abduction wedge

Fit: Adjustable straps and positioning

Clinician-directed fitting and wear required.

Product details & assistance

For fitting or recovery questions, contact your healthcare professional. For product or warranty assistance, contact Dr. Franklyn's with SKU **SS-04** and proof of purchase.

SCAN FOR PRODUCT DETAILS



drfranklyns.com

Distributed by

United Yoram Inc.

Brooklyn, NY 11236

drfranklyns.com

support@drfranklyns.com

Instagram: @drfranklyns

SS-04 | Rev. 2.0 | October 2026



SHOULDER SUPPORT / SKU SS-04

Deluxe Shoulder Immobilizer Brace

WITH ADJUSTABLE ABDUCTION WEDGE



USER MANUAL

Read all instructions before use.

Use as directed by your healthcare professional.

SUPPORT | STABILITY | COMFORT

Prepare & position the arm

1 PREPARE THE BRACE

Sit securely. Keep the injured arm supported by a helper. Open the waist, shoulder and arm straps. Have the assembly configured for the injured left or right arm.

2 POSITION THE BODY SUPPORT

Place the padded body panel against the side of the torso. Fasten the waist strap snugly. The wedge and arm support must match the position prescribed by your clinician.

3 SUPPORT THE FOREARM

With help, rest the forearm on the support platform and seat the elbow securely. Keep the hand and wrist in the position directed by your clinician.

Do not lift or rotate the injured arm beyond your prescribed limits during fitting.

Secure & check the fit

4 FASTEN THE ARM STRAPS

Secure the upper-arm, forearm and wrist straps evenly. The brace should feel secure without causing pressure, pain or restricted circulation.

5 ADJUST THE SHOULDER STRAP

Route the padded strap over the opposite shoulder and attach it to the brace. Adjust its length to support the arm without pulling on the neck.

6 CONFIRM THE WEDGE SETTING

Have your clinician confirm the abduction wedge position. Do not change the angle or remove the wedge without medical guidance. Recheck every fastener.

Comfort & circulation: Fingers should stay warm and their usual color. For numbness, tingling, swelling, cold fingers, discoloration or increased pain, support the arm, loosen tight straps and seek prompt medical advice.

Protect your arm. Care for your brace.

IMPORTANT SAFETY

Use under the direction of a qualified healthcare professional. Follow the prescribed arm position and wearing schedule, including instructions for sleep, bathing and exercise. Do not drive or operate machinery while immobilized.

Do not apply over open wounds, infected skin or an unprotected incision unless directed by your clinician. Seek fitting advice if you have impaired circulation or sensation. Check skin daily for irritation or pressure marks.

CLEANING & INSPECTION

Follow the care label. Close hook-and-loop fasteners. Hand wash removable fabric parts in cool water with mild soap. Do not soak rigid support parts. Air dry completely; do not bleach, iron or tumble dry. Reassemble as fitted by your clinician.

Check straps, fasteners and rigid parts regularly. Stop using damaged parts and contact product support.