

PRODUCT INFORMATION

Keep this guide for future reference. Read all instructions before use and follow your healthcare professional's directions when the product is prescribed.

CUSTOMER SUPPORT

Distributed by United Yoram Inc. | Brooklyn, NY 11236
drfranklyns.com | support@drfranklyns.com
Instagram: @drfranklyns

DOCUMENT CONTROL

Brand: Dr. Franklyn's | Model: SS-1 | Revision: 1.0 | Issue date: September 2026

SCAN FOR PRODUCT DETAILS



drfranklyns.com

INTENDED USE AND MEDICAL PURPOSE

The SS-1 is a reusable shoulder orthosis intended to support the forearm and help immobilize the shoulder by securing the arm against the torso. They may be used, when directed by a healthcare professional, following shoulder dislocation, rotator-cuff injury, sprain or strain, fracture management, or surgery. Available in gray and green variants.

SUPPORT PROVIDED

- Torso band and upper-arm cuff help restrict shoulder and upper-arm movement.
- Forearm and wrist supports maintain the arm in a supported position.
- Adjustable shoulder and body straps permit left- or right-side fitting.
- Lightweight breathable materials support clinician-directed daily wear.

DO NOT USE OR SEEK ADVICE

- Do not use without medical evaluation after severe trauma, suspected fracture or dislocation, deformity, or sudden loss of arm function.
- Do not apply over open wounds, infected skin, or an unprotected surgical incision.
- Seek fitting advice if you have impaired circulation, sensation, swelling, or a condition requiring a specific arm angle.



HEALTH | SUPPORT | RECOVERY



MODEL SS-1

UNIVERSAL SHOULDER
IMMOBILIZER

ADJUSTABLE ARM AND SHOULDER STABILIZATION

READ INSTRUCTIONS COMPLETELY BEFORE USE

BEFORE FITTING

- Confirm the correct side and inspect all cuffs, straps, stitching, and closures.
- Sit upright. Keep the injured arm supported while opening and positioning the immobilizer.

FITTING INSTRUCTIONS

- 1

TORSO
Wrap the wide body band around the torso at a comfortable level.
- 2

FOREARM
Place the forearm across the abdomen in the position directed by your clinician.
- 3

CUFFS
Secure the upper-arm and wrist/forearm cuffs without compressing nerves or circulation.

FIT CHECK

The elbow should be supported, the arm should rest close to the torso, and the shoulder strap should not press on the neck. Fingers must remain warm, normally colored, and freely movable.

COMPLETE THE FITTING

- 4

SHOULDER
Route the long strap over the opposite shoulder and attach it without neck pressure.
- 5

STABILIZE
Fasten the arm assembly to the torso band so the elbow remains close to the body.
- 6

RECHECK
Confirm the hand, wrist, elbow, and shoulder remain supported and comfortable.

USE AND WEAR

- Follow the arm position and wear schedule prescribed by your healthcare professional.
 - Remove only as permitted for skin checks, hygiene, and prescribed exercises. Support the arm during removal.
- Loosen or remove the device and seek advice for increasing pain, swelling, numbness, tingling, cold or discolored fingers, skin injury, breathing difficulty, or loss of hand movement.

IMPORTANT SAFETY INFORMATION

- Incorrect arm position may stress healing tissue. Follow postoperative or fracture precautions exactly.
- Do not drive, operate machinery, or perform activities requiring two hands while immobilized.
- Single-patient use only. Do not cut straps or alter the construction.

CLEANING AND STORAGE

- Fasten hook-and-loop closures before washing. Hand wash in cool water with mild soap.
- Do not bleach, wring, machine dry, iron, or dry clean.
- Air dry completely away from heat. Store clean and dry with straps untangled.

WARRANTY AND EXPECTED LIFE

One-year limited warranty from the original purchase date. Expected product life: three years under normal fitting, inspection, care, cleaning, and storage. Replace sooner if damaged, stretched, or unable to provide the intended support. Retain proof of purchase.