

PRODUCT INFORMATION

Keep this guide for future reference. Read all instructions and follow your healthcare professional's directions when the brace is prescribed.

CUSTOMER SUPPORT

Distributed by United Yoram Inc. | Brooklyn, NY 11236
drfranklyns.com | support@drfranklyns.com
Instagram: @drfranklyns

DOCUMENT CONTROL

Brand: Dr. Franklyn's | Model: TLSO | Revision: 1.0 | Issue date: September 2026

SCAN FOR PRODUCT DETAILS



drfranklyns.com

INTENDED USE AND MEDICAL PURPOSE

This reusable thoracic-lumbar-sacral orthosis is intended to provide external trunk support, encourage upright posture, and help limit unwanted thoracic and lumbar motion. It may support clinician-directed management of postural weakness, back discomfort, injury recovery, or postoperative care.

SUPPORT PROVIDED

- Shoulder straps provide adjustable postural guidance across the upper torso.
- The waist section supplies circumferential lumbar support and stabilization.
- Breathable, durable materials permit adjustable daily wear.
- Size variants are available in L, XL, and XXL.

DO NOT USE OR SEEK ADVICE

- Do not use over open wounds, infected skin, or an unprotected surgical incision.
- A clinician should select and fit the brace after spine surgery, fracture, or with significant spinal deformity.
- Do not overtighten the shoulder or waist straps or use if breathing, circulation, or sensation is affected.



HEALTH | SUPPORT | RECOVERY



MODEL TLSO

POSTURE CONTROL
TLSO BRACE

THORACIC AND LUMBAR SPINE SUPPORT

READ INSTRUCTIONS COMPLETELY BEFORE USE

BEFORE FITTING

- Inspect the brace and straps. Wear a thin shirt beneath the brace to protect the skin.
- Loosen the waist and shoulder straps before positioning.

FITTING INSTRUCTIONS

- 1

POSITION
Center the rear support along the spine with the lower section across the lumbar area.
- 2

SHOULDERS
Place each shoulder strap over the corresponding shoulder without twisting.
- 3

WAIST
Wrap the waist belt evenly and fasten it at the front while standing upright.

FIT CHECK

The rear panel should follow the center of the back, the waist support should remain level, and the shoulder straps should guide rather than force posture.

COMPLETE THE FITTING

- 4

GUIDE
Draw the shoulder straps down and back to encourage posture without pulling the shoulders painfully.
- 5

BALANCE
Adjust both sides evenly so the support remains centered and level.
- 6

RECHECK
Breathe, sit, and move gently. Loosen any strap that causes pressure or arm symptoms.

USE AND WEAR

- Use only for the duration and activities directed by your healthcare professional.
- Remove periodically to inspect the shoulders, underarms, waist, and back for pressure areas.

Remove the brace and seek medical advice for increasing pain, numbness, tingling, unusual swelling, skin breakdown, or pale, blue, cold, or discolored skin.

IMPORTANT SAFETY INFORMATION

- Never use the shoulder straps to force the torso beyond a comfortable upright position.
- Seek advice for new weakness, radiating pain, numbness, breathing difficulty, or loss of balance.
- Single-patient use only. Do not cut, bend, or modify structural components.

CLEANING AND STORAGE

- Fasten all hook-and-loop closures before cleaning. Remove detachable rigid stays only when the care label permits.
- Hand wash in cool water with mild soap. Do not bleach, wring, machine dry, iron, or dry clean.
- Air dry completely away from direct heat and sunlight. Store flat or loosely folded in a clean, dry place.

WARRANTY AND EXPECTED LIFE

One-year limited warranty from the original purchase date. Expected product life: three years under normal fitting, inspection, care, cleaning, and storage. Replace sooner if damaged, stretched, or unable to provide the intended support. Retain proof of purchase.