



CUSTOMER SUPPORT

Distributed by United Yoram Inc. | Brooklyn, NY 11236
drfranklyns.com
support@drfranklyns.com
Instagram: @drfranklyns

PRODUCT INFORMATION

Product: Water Circulation Therapy
Model / SKU: WCP-2
User manual

WHAT IS INCLUDED

- 1 Cooler
- 1 Pump
- 1 Therapy pad
- 1 Connecting hose

KEEP THIS GUIDE

Keep this manual for future reference. Read all instructions and safety warnings before use. Follow your physician's directions for treatment.

IMPORTANT SAFETY INFORMATION

Check the treated skin every 5-10 minutes during therapy. Stop immediately for pale or bluish skin (cyanosis), excessive redness, burning, severe stabbing pain or numbness.

DO NOT USE

- Severe peripheral vascular disease or regional blood circulation disorders.
- For cold use: Raynaud's disease, cold allergy or cold urticaria.
- Loss of local skin sensation or severe numbness.
- Large open wounds or infectious ulcers at the treatment site.

PROTECT SKIN

Place clean non-woven fabric, gauze or thin clothing between the pad and skin. Never allow the water pad to touch bare skin, to help prevent cold injury or burns. For warm therapy, use comfortably warm water; never use boiling or scalding water.

HEIGHT LIMITATION

Do not place the bucket more than 38 cm (15 inches) above the pad. Excess height can increase water pressure in the pad and cause discomfort.

ELECTRICAL & PUMP SAFETY

- Do not start the pump without water in the bucket. Dry running may overheat and damage the motor.
- Do not immerse the lid power port or power cord in water.
- Never plug in or unplug the adapter with wet hands.

STOP USE

Stop treatment immediately if any of the skin changes or symptoms listed above occur. Switch off and unplug the system.



WATER CIRCULATION THERAPY

HOT & COLD WATER CIRCULATION



MODEL WCP-2

USER MANUAL

READ ALL INSTRUCTIONS BEFORE USE

HEALTH | SUPPORT | RECOVERY

INTENDED USE & BENEFITS

The system circulates water through a therapy pad for hot or cold application to the target area. Cold therapy may help relieve pain and swelling. Warm therapy may help ease muscle aches and stiffness.

WHEN IT IS USED

Cold: May be used following acute soft tissue injuries, strains, sprains or orthopedic procedures, as directed by your clinician.

Warm: May be used for muscle discomfort and stiffness. Do not apply heat to a new injury or an area that is swollen, bruised or has an open wound.

HOW THE SYSTEM WORKS

The pump circulates the water added to the bucket through the connecting hose and therapy pad. Secure the pad comfortably with its straps. Select warm or cold water for the intended therapy.

OPERATING REFERENCE

Cold use: ice-to-water	Approx. 1:1
Cold preparation	3-5 minutes
Cold session	15-30 minutes*
Warm session	As directed by clinician
Skin checks	Every 5-10 minutes
Max. bucket height above pad	38 cm / 15 inches

*Or follow your physician's instructions.

SETUP AND OPERATION

1. Prepare the water
Cold use: Open the lid. Add cold water to the marked Water Line, then add ice (approximately 1:1 ice to water). Secure the lid. Let stand for 3-5 minutes, then gently shake to lower the water temperature.

Warm use: Add comfortably warm water to the marked Water Line. Do not add ice. Never use boiling or scalding water. Secure the lid.

2. Apply the therapy pad
Place clean non-woven fabric, gauze or thin clothing over the affected area. Do not let the water pad contact bare skin.

Wrap the unfilled pad smoothly around the target area, such as the knee or shoulder. Secure with hook-and-loop straps. The fit should be snug, but not too tight.

3. Connect and start
Connect the delivery hose between the bucket outlet and the pad port. A click confirms each quick-connect fitting is fully locked.

Plug the power adapter into a wall outlet and connect its other end to the power port on the bucket lid. Turn on the switch. The pump circulates water through the tubing and therapy pad.

4. End therapy
Follow your clinician's directions for the therapy type, water temperature and session duration. For cold use, the recommended session is 15-30 minutes, or as directed by your physician. Turn off and unplug after treatment.

MONITOR DURING TREATMENT

Check skin every 5-10 minutes and follow all safety information on the reverse side.

AFTER EACH USE

- Switch off and unplug the system.
- Press the quick-release button to disconnect the hose.
- Open the drain port on the therapy pad to empty residual water.
- Pour out the remaining water from the bucket.

CLEANING & MAINTENANCE

Rinse the inside of the bucket and connecting hose with warm water mixed with neutral detergent. Flush thoroughly with clean water.

DRY THOROUGHLY

Open the bucket lid. Place the therapy pad in a cool, well-ventilated area and allow it to air dry naturally to help prevent mold growth inside.

STORAGE

Store in a cool, dry place. Avoid direct sunlight and heavy pressure.

ALWAYS ADD WATER FIRST

Never power on or start the pump with an empty bucket. Dry running can overheat and damage the motor.